

388 XL

Mid size scooters



Your new route to independence

RASCAL 

388 XL

388 XL Specification:

Type Class	C / 3
Overall length	125cm / 49"
Overall width	54cm / 21"
Seat width	45.7cm / 18"
Seat depth	45.7cm / 18"
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 33Ah x 2
Maximum speed	9.65km/h / 6mph
Ground clearance	13cm / 5"
Maximum kerb height	5cm / 2"
Range*	up to 32km / 20 miles
The maximum safe slope	8°
Turn-around width	1.7m / 5'6"
Colour	Black / Red

*Range on full charge and flat ground

Rascal 388 XL Profile:

The Rascal 388 XL is the perfect Class 3 scooter with 6 mph top speed, stylish 11" wheels, full road lighting and with its narrow chassis it combines safety with performance and manoeuvrability. Sturdier in build and ideal for people requiring additional comfort on those slightly longer journeys - a superb scooter for days out and leisure trips. With a comfortable high back seat with headrest, width adjustable armrests and tiller, this scooter can be 'adapted' to suit the individual's requirements. 388 S and 388 Deluxe models also available.

Local Rascal Mobility Retailer:

Main Features:

- Outstanding reliability
- Flip-up width adjustable armrests
- Comfortable high back, height adjustable, sliding swivel seat
- Narrow compact chassis
- Lift-up tiller head for easy access
- Captain's seat
- Smooth controls
- One wing mirror as standard
- Front basket as standard

DISCOVER

Other Models:

- Enquire about the Rascal 388 S and 388 D today - narrow 4mph scooters for the everyday.



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E - v 6 May 2019