

Rascal 388 Range



Pictured below
388 XL in black



Optional
canopy
available

Rascal 388 Range



Rascal 388 Standard
available in red and silver



Rascal 388 Deluxe
only available in red



Rascal 388 XL
available in black and red

Rascal 388 Range Specification:

	Standard	Deluxe	XL
Overall length	127cm / 50"	127cm / 50"	127cm / 50"
Overall width	56cm / 22"	56cm / 22"	56cm / 22"
Batteries	2x34amp	2x34amp	2x34amp
Ground Clearance	10cm / 4"	10cm / 4"	13cm / 5"
Maximum range*	32km / 20 miles	32km / 20 miles	32km / 20 miles
Maximum speed	6.4kph / 4mph	6.4kph / 4mph	9.8kph / 6mph
Maximum rider weight**	136kg / 21 stone	136kg / 21 stone	136kg / 21 stone
Total weight	74kg / 11 stone	85 kg / 13 stone	85 kg / 13 stone

*Range is calculated under test conditions.
Weight of rider, terrain and weather conditions can affect range of travel.
**Maximum rider weight is reduced to 91kg (200lbs) when climbing a 15% incline.

Independence is priceless. That is why every scooter in the Rascal 388 range offers you the choice and freedom you deserve. Your daily life will take on a whole new perspective as many more destinations become possible with a Rascal 388S, 388D or 388XL. When it comes to choosing a mid-size scooter with proven performance, your expectations will be exceeded!

With three models to choose from, you can own the ideal scooter for your needs. All are designed with a narrow, compact chassis for ease of manoeuvring and a sliding, swivel seat for comfort and accessibility.



All seats are
sliding with swivel



Carrying capacity
of up to 300lbs



Narrow compact
chassis



www.facebook.com/RascalMobility
www.twitter.com/RascalMobility

Electric Mobility Euro Ltd, Canal Way, Ilminster, Somerset, TA19 9DL
0800 252614 | sales@electricmobility.co.uk | www.electricmobility.co.uk

Due to a policy of continual improvements, Electric Mobility reserve the right to change specification without prior notice.

Local Rascal Mobility Retailer: