

# Ultralite 480




Voted

**BEST  
BUY**

in Which? Magazine



Lightweight  
Components



Good  
Legroom

# Ultralite 480 Voted ‘Best Buy’ in the Mobility Scooter Reviews by Which? Magazine

## Ultralite 480 Specification:

Type Class	A / 2
Overall length	102cm / 40"
Overall width	51cm / 20"
Floor to seat height	45-55cm (min-max)
Seat width	43.1cm / 17"
Seat depth	33cm / 13"
Maximum carrying capacity	113kgs / 18 stone
Battery voltage and capacity	12V / 12Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	5cm / 2"
Maximum kerb height	4cm / 1.5"
Range*	up to 13km / 8 miles
The maximum safe slope	6°
Turn-around width	1.62m / 5'4"
Colour	Black / Red / Teal

\*Range on full charge and flat ground.

## Main Features:

- Car transportable
- Fold down backrest
- Flip-up width adjustable armrests
- Removable front basket
- Tiller adjustment
- Solid tyres all-round
- Charger storage in battery box
- Swivel seat
- Single heaviest component (front chassis) 13kgs / 29lbs

## Optional Extras:

- A cane and crutch holder bag



## Ultralite 480 Profile:

The Ultralite 480 is the ultimate transportable scooter. With its compact and lightweight design, the 480 ensures you stay independent at all times. This scooter has a proven reputation for outstanding design and reliability. Easy to transport, the 480 can be dismantled in seconds; no cables to worry about, just quick-release locks, convenient lifting handles and lightweight manageable components. Easy to stow, the 480 can be stored vertically in your home if space is an issue, and fits easily into most car boots when taken apart.



## What the Which? review states...

“ We were so impressed with this mobility scooter that we’ve made it a Best Buy. It gives a reassuringly stable ride, manoeuvres well around corners, and is easy to get on and off. You won’t have to recharge the battery too often, and you shouldn’t have too much trouble lifting its component parts into the boot of a car when it’s dismantled. ”

taken from Which? magazine review

## Which? Ratings

Portability	★★★★☆
Ease of handling	★★★★☆
Control and manoeuvrability	★★★★☆
Getting on and off	★★★★★
Using controls	★★★★★

Score  
79%



## What’s it great at?

You’ll get a comfortable ride on smooth surfaces, and although the tiller vibrates when driven over gravel the padded seat ensures the driver doesn’t get too shaken up. It slows down considerably when going uphill but good stability means you’ll always feel secure.

The wheels are very responsive so it’s very easy to manoeuvre around obstacles, and we like the simple control panel with its clear and intuitive illustrations. The battery life indicator is clearly visible, so you’ll have plenty of warning before you run out of power.

It’s easy to assemble, and a loud click will let you know when the parts are fitted into place correctly.

## Can I rely on the battery to go the distance stated on full charge?

It takes less than six and a half hours to fully charge the battery, which is a few hours quicker than some we’ve tested. It can be charged separately or when it’s in place, and is easy to connect to the charger.

## Should I buy it?

Yes, this Best Buy scooter is easy to steer, easy to fit into a car boot and gives a stable, comfortable ride.

# Ultralite 480

**With 3 fabulous colours to choose from!  
What colour will you buy?**

★ **Black**



★ **Red**



★ **Teal**



Local Rascal Mobility Retailer:



Specifications subject to change without prior notice  
E. and O.E- v 4 Final-August 2014