

# 388 S & 388 Deluxe

Mid size scooters



Your new route to independence

**RASCAL** 

# 388 S & 388 Deluxe

## 388 S & D Specification:

Type Class	B / 2
Overall length	125cm / 49"
Overall width	54cm / 21"
Floor to seat height (cm)	56-64 <b>(S)</b> / 57-65 <b>(D)</b>
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16" <b>(S)</b> / 45.7cm / 18" <b>(D)</b>
Maximum carrying capacity	136kgs / 21 stone
Battery voltage and capacity	12V / 33Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	10cm / 4"
Maximum kerb height	5cm / 2"
Range* - standard battery pack	up to 32km / 20 miles
The maximum safe slope	8°
Turn-around width	1.7m / 5'6"
Colour - 388 S	Red / Silver
Colour - 388 Deluxe	Red

\*Range on full charge and flat ground

## Rascal 388 S & D Profile:

The Rascal 388 S and 388 Deluxe are both well established, great value for money and suit your everyday mobility needs. At 21" wide, it is the narrowest width of any scooter in its class & is equally at home inside a busy supermarket as it is on the pavement. With a comfortable high back seat with headrest, smooth controls and adjustable tiller, these high performance, reliable scooters are your ideal mobility partner.

Local Rascal Mobility Retailer:

## Main Features 388 S:

- Outstanding reliability
- Flip-up width adjustable armrests
- Comfortable high back, height adjustable, sliding swivel seat
- Narrow compact chassis

## Main Features 388 D:

- All of the above features of the 388 S
- Lift-up tiller head for easy access
- Captain's seat
- Smooth controls



UPGRADE

## Upgrade:

- The Rascal 388 XL, offers all the great features of the 388 Deluxe, plus a maximum speed of 6 mph and a ground clearance of 13cm



Specifications subject to change without prior notice  
Images are for illustrative purposes only  
E. and O.E - v 5 February 2019