



VP172AA

# **Lightweight Tri Walker with Bag and Basket**

The aluminium Tri-Walker from Aidapt is sturdy yet easy to manoeuvre, featuring loop brakes so the tri-walker can be stopped (by squeezing) and locked (by pushing down). Suitable for both indoor and outdoor use, the Aidapt Tri-Walker with Bag and Basket folds away for transport and storage. Featuring a unique cable guide to eliminate the need for plastic cable ties, with height adjustable handles for comfort, the three-wheeled walker is also supplied complete with a wipe-clean, specially shaped tri-walker bag and basket with tray.

- *Ideal for indoor and outdoor use*
- *Adjustable height handles for a comfortable fit*
- *Easy folding for storage and transport*
- *Simple-to-use lock brakes and cable tidy for safety*
- *Supplied with bag, basket and tray*

Height (mm)	830-990 (7 height adjustments)
Net weight (with bag, basket & lid)	6.8
Maximum User Weight (kg)	132
Maximum User Weight (stones)	21