

The  
**Jubilee**

Tilt-in-Space Dual-motor

Made  
in the UK









# Jubilee

## Lateral Back

Choose our furniture with confidence -  
All riser recliner chairs are tested by FIRA -  
The leading authority in furniture safety by  
design.



### Riser Recliners available from STOCK in the following fabrics



2 or 3-Seater static  
sofa available  
Made to order



Button & Mini Lateral Back  
Made to order or Made to measure

Colours may vary slightly from images shown which are for illustrative purposes only.

The Jubilee is our UK Premium Tilt-in-Space (TIS) Dual-motor Cosi Chair designed to maximise support throughout motion with a range of stylish back designs & upholstery options to provide exceptional postural seating positions to suit all individuals. The Jubilee is available in Waterfall or Lateral Back from stock with other back options & fabrics available Made to order.

### Main Features:

- Premium Seminar TIS mechanism with Variable Angle Lift\* (VAL\* - M & L only. Ask for more details)
- Arm covers and Antimacassar included
- 5-button power on handset
- Choice of configurable back cushions
- Heat & massage available Made to order
- Drop in seat cushions for pressure care issues available Made to order. (Ask for more details.)
- 2 or 3-Seater static sofa available Made to order
- 5 Year Guarantee on frames, mechanisms and all electrics
- Quality Craftsmanship, **BRITISH MADE**

Local Cosi Chair Retailer:

### Jubilee Lateral Size Guide:

\* Also available in Petite - measurements Seat H 16" x Seat W 18" x Seat D 16" - on Made to order

	<b>Small</b>	<b>Medium</b>	<b>Large</b>
Seat height	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Seat width	45.7cm / 18"	50.8cm / 20"	55.9cm / 22"
Seat depth	45.7cm / 18"	50.8cm / 20"	53.3cm / 21"
Arm height	19.1cm / 7½"	19.1cm / 7½"	19.1cm / 7½"
Back height	68.6cm / 27"	73.7cm / 29"	78.7cm / 31"
Overall width	76.2cm / 30"	81.3cm / 32"	86.4cm / 34"
Overall height	108cm / 42½"	110.5cm / 43½"	113cm / 44½"
Weight limit	20 stone / 127kg	25 stone / 158.8kg	25 stone / 158.8kg

