

## Technical Specification - Rascal 320 Compact

Model No.	P320 Specifications
Length	32.5" (825mm)
Width	23" (580mm)
Seat Width	18" (460mm)
Seat Height(from deck)	15"~17" (380~430mm)
Seat Height(from ground)	18"~20" (460~510mm)
Speed	6 kph (3.73 mph)
Weight Capacity	113.5 kg (250 lbs)
Motor	DC24V, 120W, 4000rpm
Controller	Dynamic Shark
Battery	12V, 20Ah * 2 pcs
Charger	3A off board
Gradient	8°
Front Wheel	7*1-3/4" Pneumatic or PU
Rear Wheel	9" Pneumatic or PU
Total Weight (with battery)	59.6 kg (131 lbs)
Battery Weight	14.2 kg (31 lbs)
Seat Weight	18.2 kg (40 lbs)
Chassis Weight	27.2 kg (60 lbs)