

Rascal P312 Turnabout SL

Mid-range powerchairs



Your new route to independence

RASCAL 

P312 Turnabout SL

P312 Turnabout SL Specification:

Type Class	B / 2
Overall length	96cm / 38"
Overall width	61cm / 24"
Floor to seat height	57- 69.5cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm / 16"
Maximum carrying capacity	159kgs / 25 stone
Battery voltage and capacity	12V / 50Ah x 2
Maximum speed	6km/h / 3.8mph
Ground clearance	5cm / 2"
Maximum kerb height	5cm / 2"
Range* - standard battery pack	up to 29km / 18 miles
The maximum safe slope	10°
Turn-around width	1.18m / 3'10"
Colour	Red

*Range on full charge and flat ground

Main Features:

- 5 inch Column Seat Lift
- Front-Wheel or Rear-Wheel Drive
- Flip-up width adjustable armrests
- Convenient joystick charging port
- Swing-away joystick
- High back adjustable seat
- Fold-away footplates front and rear

**OPTIONAL
ACCESSORY**



Optional Extras:

- Option of swing-away legrests
- Option of elevating legrests

Rascal P312 Turnabout SL Profile:

The Rascal P312 Turnabout with Seat Lift gives flexibility to the user including a dual directional drive with a rotating seat giving you the option of Front-Wheel Drive for superb manoeuvrability indoors or Rear-Wheel Drive for extra traction outdoors. The P312 Turnabout Seat Lift offers well defined driving positions with a swing-away joystick and a fully adjustable seat and headrest for day-long comfort.

Local Rascal Mobility Retailer:



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E. - v 5 Final-May 2019