

Pioneer

Performance scooters



Your new route to independence

RASCAL 

Pioneer

Pioneer Specification:

Type Class	C / 3
Overall length	147cm / 58"
Overall width	64.5cm / 25"
Floor to seat height	65-68cm (min-max)
Seat width	45.7cm / 18"
Seat depth	45.7cm / 18"
Maximum carrying capacity	159kgs / 25 stone
Battery voltage and capacity**	12V / 75Ah x 2
Maximum speed	12.9km/h / 8mph
Ground clearance	9cm / 3.5"
Maximum kerb height	10cm / 4"
Range* - standard battery pack)	up to 56km / 35 miles
The maximum safe slope	10°
Turn-around width	2.02m / 6'7"
Colour	Black / Red

*Range on full charge and flat ground

Rascal Pioneer Profile:

The Rascal Pioneer is a proven performer with a fantastic reliability and durability record. With full coil suspension and a 35 mile range its powerful motor makes it an ideal scooter for the more challenging environments. It features a delta-bar tiller with adjustment, a sliding swivel seat with headrest and padded arms as well as showerproof switches and a large basket making the Pioneer one of the most popular scooters on the market.

Local Rascal Mobility Retailer:

Main Features:

- Front and rear coil suspension
- Removable front basket
- Gas strut tiller adjustment with useful battery charging port
- Modern delta handlebars
- Showerproof panel and twin mirrors
- Large shock absorbing tyres for maximum comfort on the move
- Sliding swivel seat with armrests and headrest for maximum comfort

OPTIONAL
ACCESSORY

Optional Extras:

- A stylish rear mounted lockable box fitted on a sliding mount to enable easy access without adjusting the seat
- An exclusive Rascal basket bag
- Standard all weather canopy - upgrade to a polycarbonate front shield for an additional cost



Specifications subject to change without prior notice
Images are for illustrative purposes only
E. and O.E - v 8 March 2019