



Your new route to independence



Veo Sport SR

Veo Sport SR Specification:

•	
Type Class	B/2
Overall length	108cm / 42.5"
Overall width	60cm / 23.6"
Floor to seat height	47-49cm (min-max)
Seat width	45.7cm / 18"
Seat depth	40.6cm/16"
Maximum carrying capacity	130kgs / 20.5 stone
Battery voltage and capacity	12V / 22Ah x 2
Maximum speed	6.4km/h / 4mph
Ground clearance	6cm / 2.4"
Maximum kerb height	5cm / 2"
Range* - 22Ah battery pack	up to 24km / 14.9 miles
The maximum safe slope	6°
Turn-around width	1.40m / 4'6"
Colour	Neon Blue / Mulberry
Single Heaviest component 20.2kgs / 44.5lbs	

Main Features:

- Car transportable
- All-round suspension
- Stylish sports seat
- Tiller charging point
- LED front light
- Front bumper
- Dual compound puncture proof tyres
- Comfortable delta tiller
- Rear safety reflectors

OPTIONAL



NEW For 20

- Available with lightweight 14Ah battery pack (Range of up to 15.3km / 9.5 miles)
- Rascal arm rest bag



Rascal Veo Sport SR Profile:

*Range on full charge and flat ground

MOBILIT

It's all about the drive. The Rascal Veo Sport SR is fully loaded with features in 2 striking colours. A modern, sporty design with class leading all-round suspension. With 22Ah batteries, a superb capacity for a scooter of its size, affording a range of up to 14.9 miles. A padded delta tiller boasting a convenient charging point and luxury seat, detailed with unique white accent stitching, the Veo Sport SR is a real pleasure to ride.

The Veo Sport SR's easy and robust take apart mechanism allows you to store or drive away with it in the boot of your car, so it is very practical yet with its distinctive looks, exclusive all-round suspension plus LED front light, this is the travel scooter with the edge.



